

Example Report

Tumor Monitoring Blood Test Report
腫瘤監測血液測試報告

Referral: 轉介:	iCare Medical Centre	Reference No.: 參考編號:	---
Examinee Name: 受測者姓名:		Specimen No.: 樣本編號:	
ID No./Passport No.: 身份證明文件號碼:	---	Specimen Type: 樣本類別:	Blood
D.O.B: 出生日期:		Specimen Collection Date: 樣本抽取日期:	26 June 2019
Sex: 性別:	Male	Specimen Received Date: 樣本接收日期:	28 June 2019
Race: 種族:	Caucasian	Reporting Date: 報告日期:	3 July 2019
Clinical History/Referral Reason: 臨床病歷 / 轉介原因:	PSA 1.3, History of Prostate cancer about 10 years ago	Report No.: 報告編號:	CTM-35IR6N

Test Items 檢驗項目	Results 結果
<u>Circulating Tumor Cell (CTC) Counting*</u> <u>血液循環腫瘤細胞檢測*</u>	Favorable 陰性 2 CTCs / 7.5mL of Blood 血液循環腫瘤細胞 / 7.5毫升血液

CTCs were identified by immuofluorescence staining. 利用免疫熒光染色法檢測血液循環腫瘤細胞。
 Cytokeratin(+) EpCAM(+) and CD45(-) cells are defined as CTCs. Cytokeratin(+) EpCAM(+) 及 CD45(-) 細胞被定義為血液循環腫瘤細胞。
 The identified CTCs are shown in Appendix of the report. 檢測出來的血液循環腫瘤細胞顯示在本報告附錄上。

MT102523

Approved Signatory: 
 Registered Medical Laboratory Technologist (Part I)

* Unfavorable result if CTC number is ≥ 10 . Favorable result if CTC number is ≤ 9 .
 * 血液循環腫瘤細胞檢測數量 ≥ 10 為陽性結果， ≤ 9 為陰性結果。
 Note 1: This report is intended for information only and serves no purpose for diagnosis and/or treatment.
 Note 2: This report must always be interpreted by physician in connection with correlated clinical conditions of the examinee.
 Note 3: The sensitivity and specificity of the above CTC test is 72.2% and 93.7%, respectively, according to current study for epithelial carcinoma only.
 Note 4: This report should not be reproduced, except in full, without the written approval of Acecgt Diagnostic Limited.

Example Report


Examinee Name:
受測者姓名:
Specimen No.:
樣本編號:

Reporting Date: 3 July 2019
報告日期:
Report No.: CTM-351R6N
報告編號:

Appendix 附錄

CTCs Images 血液循環腫瘤細胞影像

The identified CTCs are shown at the center of each image below:
以下每幅影像的中心位置顯示檢測出來的血液循環腫瘤細胞:



-End of the report-